

# StepMate

## Stepper

Positive-displacement, continuously adjustable, repetitive dispenser



## Features

- Slim, light and robust
- Wide dispensing volume from 10 $\mu$ L to 5mL
- Up to 48 dispensing steps in succession
- Equipped with durable dispenser tip stop lever
- Fits different sizes of disposable polypropylene dispenser tips from 0.5mL to 50mL

## Applications

Ideal for any application requiring repetitive dispensing in the clinical, biological, chemical, food and beverage, pharmaceutical industry.

### How to select dispensing volume

Find the required dispensing volume in the table.

Use the volume selection dial to set the required volume.

Select and insert a suitable dispenser tip. Follow the maximum dispensing steps as stated.

Dispenser tip Nominal Volume

Dispensing Volume per Step

| Tips \ Sel. | 1            | 2            | 3            | 4            | 5            |
|-------------|--------------|--------------|--------------|--------------|--------------|
| 0.5mL       | 10 $\mu$ L   | 20 $\mu$ L   | 30 $\mu$ L   | 40 $\mu$ L   | 50 $\mu$ L   |
| 1.25mL      | 25 $\mu$ L   | 50 $\mu$ L   | 75 $\mu$ L   | 100 $\mu$ L  | 125 $\mu$ L  |
| 2.50mL      | 50 $\mu$ L   | 100 $\mu$ L  | 150 $\mu$ L  | 200 $\mu$ L  | 250 $\mu$ L  |
| 5.00mL      | 100 $\mu$ L  | 200 $\mu$ L  | 300 $\mu$ L  | 400 $\mu$ L  | 500 $\mu$ L  |
| 12.50mL     | 250 $\mu$ L  | 500 $\mu$ L  | 750 $\mu$ L  | 1000 $\mu$ L | 1250 $\mu$ L |
| 25.00mL     | 500 $\mu$ L  | 1000 $\mu$ L | 1500 $\mu$ L | 2000 $\mu$ L | 2500 $\mu$ L |
| 50.00mL     | 1000 $\mu$ L | 2000 $\mu$ L | 3000 $\mu$ L | 4000 $\mu$ L | 5000 $\mu$ L |
| Steps       | 48           | 23           | 15           | 11           | 8            |

## Precision Dispenser Tip

DLAB dispenser tips are convenient for viscous or foaming solutions, volatile liquids and applications where repetitive dispensing of volumes in long series is required

The dispenser tips are made from high quality materials without soft lubricants or other additives.

## Specifications

Note: The first and last dispensing liquid should be discarded for precision.

| Dial setting    |                      | 1                          | 2    | 3    | 4    | 5    |                  |              |                                                                            |
|-----------------|----------------------|----------------------------|------|------|------|------|------------------|--------------|----------------------------------------------------------------------------|
| Number of steps |                      | 48                         | 23   | 15   | 11   | 8    |                  |              |                                                                            |
| Cat. No.        | dispenser tip volume | Dispensing volume per step |      |      |      |      | Systematic Error | Random Error | Description                                                                |
|                 | mL                   | μL                         |      |      |      |      | %                | %            |                                                                            |
| 17900032        | 0.50                 | 10                         | 20   | 30   | 40   | 50   | ± 0.8            | ≤ 0.7        | 0.5mL, piston PE-HD, cylinder PP, 1 packing, 5 pcs/packing                 |
| 17900033        | 1.25                 | 25                         | 50   | 75   | 100  | 125  | ± 0.8            | ≤ 0.5        | 1.25mL, piston PE-HD, cylinder PP, 1 packing, 5 pcs/packing                |
| 17900034        | 2.50                 | 50                         | 100  | 150  | 200  | 250  | ± 0.8            | ≤ 0.4        | 2.5mL, piston PE-HD, cylinder PP, 1 packing, 5 pcs/packing                 |
| 17900035        | 5.00                 | 100                        | 200  | 300  | 400  | 500  | ± 0.4            | ≤ 0.3        | 5mL, piston PE-HD, cylinder PP, 1 packing, 5 pcs/packing                   |
| 17900036        | 12.50                | 250                        | 500  | 750  | 1000 | 1250 | ± 0.3            | ≤ 0.3        | 12.5mL, piston PE-HD, cylinder PP, 1 packing, 5 pcs/packing                |
| 17900037        | 25.00                | 500                        | 1000 | 1500 | 2000 | 2500 | ± 0.2            | ≤ 0.3        | 25mL, piston PE-HD, cylinder PP, 1 packing, 2pcs/packing with 1 pc adapter |
| 17900038        | 50.00                | 1000                       | 2000 | 3000 | 4000 | 5000 | ± 0.2            | ≤ 0.2        | 50mL, piston PE-HD, cylinder PP, 1 packing, 2pcs/packing with 1 pc adapter |

