

RE500-E&RE1000-E& RE2000-E&RE5000-E

Pilot Scale Digital Control Rotary Evaporator

Features

- LCD digital display with intuitive interface for operation and control.
- PID controller, precise control, fast setting.
- Stepless automatic motor lifting of the heat bath, smooth and silent.
- The reciprocation rotation function is the most useful for concentration of sample with power or solid particles.
- Double-layer high-efficiency condenser with anti-backflow design, large condensation surface area, strong evaporation capacity is a guarantee of high efficient sample recovery.
- Professional air tightness system design with PTFE sealing ring, featured with heat/corrosion resistance, and high durability.
- Optional vacuum controller to control the vacuum degree freely.
- Keep the outer surface of the heat bath at mild temperature with proprietary heat insulation design to avoid accidental personnel burn injury.
- Equipped with automatic switching valve, continuous collection can be carried out without affecting the negative pressure index of the system and solvent distillation.
- Optional water bath/oil bath heating mode (standard water bath).
- Solution can be continuously fed into evaporating flask.
- Digital display steam temperature sensing device can feed back the solvent steam temperature in real time



Specification

Photos				
Model	RE500-E	RE1000-E	RE2000-E	RE5000-E
Evaporating flask volume	5L	10L	20L	50L
Speed range	20-180	20-130	20-130	20-130
Display	LCD			
Evaporation Capacity (Max.L/h) For 75% Alcohol	2	3.8	5	11
Ultimate vacuum	9 mbar			
Timer	Yes			
Temperature range	Water Bath: 20- 100°C Oil Bath: 20-180°C (Optional)			
Temperature control accuracy I°C	±1			
Heating Power (W)	2000	3000	4000	5000
Temperature Display	LCD			
Lifting mode	Motorized			
Lifting Stroke (mm)	130	200	180	253
Water Bath Volume	13L	28L	40L	60L
Total Power	2200W	3200W	4200W	5200W
Time Setting Range	0-999min			
Voltage	220V			
Frequency (Hz)	50/60			
Motor	Brushless DC motor			
Heat Bath	Teflon coating			
Protection Class	IP20			
Dimension (D×W×H)	600*290*720	790*480*1130	900*490*1320	1030*510*1370
Permissible ambient temperature I°C	5-40			
Permissible ambient humidity	≤80%			
Safety features	Motor over current protection, earth leakage protection, dry heating protection, overheat protection function			

Vertical

Cat. No.	Glassware type & Volume	Cooling surface
18900674	5L Main condenser	0.28m ²
18900639	10L Main condenser	0.49m ²
18900640	10L Secondary condenser	
18900688	20L Main condenser	1.29m ²
18900689	20L Secondary condenser	
18900648	50L Main condenser	1.75m ²
18900649	50L Secondary condenser	



Evaporating flask

Cat. No.	Volume	Size
18900677	5L	RE500-E
18900643	10L	RE1000-E
18900690	20L	RE2000-E
18900650	50L	RE5000-E



Receiving flask

Cat. No.	Volume	Size
18900678	3L	RE500-E
18900644	5L	RE1000-E
18900691	10L	RE2000-E
18900651	20L	RE5000-E



Three way bottle

Cat. No.	Volume	Size
18900676	—	RE500-E
18900642	—	RE1000-E/RE2000-E/RE5000-E



Rotation axis

Cat. No.	Volume	Size
18900679	—	RE500-E
18900645	—	RE1000-E/RE2000-E/RE5000-E



Valve accessories

Cat. No.	Description	Size	Photograph
18900675	Feeding valve	RE500-E	
18900641	Feeding valve	RE1000-E, RE2000-E, RE5000-E	
18900675	Discharge valve	RE500-E	
18900641		RE1000-E, RE2000-E, RE5000-E	
18900680	Non pass valve	RE500-E	
18900681	Mixing valve	RE500-E	
18900647	Vent valve	RE1000-E, RE2000-E, RE5000-E	

Other accessories

Cat. No.	Description	Size	Photograph
18900683	5L Locking clamp 10L/20L/50L Locking clamp	30#	
18900682		40#	
18900655		50#	
18900656		60#	
18900657		80#	
18204035	5L Vacuum sealing ring	Teflon + fluorine rubber	
18203979	10L/20L/50L Vacuum sealing ring	RE1000-E, RE2000-E, RE5000-E	
18203980	Vacuum seal assembly for evaporating flask	RE1000-E, RE2000-E, RE5000-E	
18204228	Flask holder (50mL-1L and 1L-20L)	90mm	
		160mm	
18901003	Evaporation bottle bracket	adapted to RE500-E	
18901004		adapted to RE1000-E	
18901005		adapted to RE2000-E	
18901006		adapted to RE5000-E	